


THE GLOBAL
A-SERIES™
DENTAL MICROSCOPE



GLOBAL
SURGICAL™ CORPORATION
We make it easier



AXIS™ CONTROL SYSTEM

- Optional Carr Extender 1
- Inclinable Binoculars 2
- Eyecups 3
- Eyepieces 4
- Diopter Adjustment 5
- Magnification Changer 6
- Fine Focus Adjustment 7
- Optional Multi-Focal Lens 8

ALL NEW DESIGN	Sleek, innovative style for the modern dental office
INTUITIVE CONTROL	The new AXIS™ Control System makes it easy to position and control
BRIGHTER LIGHT	With over 100,000 Lux output, we offer the brightest LED light source available
ENHANCED COMFORT	Updated ergonomic design to promote long-term health benefits

RESTORATIVE DENTISTRY | ENDODONTICS

THE MICROSCOPE EXPERIENCE

Global is committed to delivering the best microscope experience available. The A-Series™ microscope is designed for intuitive control while significantly improving visualization for earlier diagnosis and treatment, setting a new standard in dental microscopy.

- Available with 3, 4, or 6-steps of magnification
- AXIS™ Control System delivers magnification changer, ergonomic handles, and tilt tension controls at one location bilaterally
- Updated binoculars with diopter eyepieces
- Brightest LED light source available
- Three-stage filter changer: clear, amber, and green
- Custom mounting options to suit any office space: floor, wall or ceiling
- Limited lifetime warranty
- The only dental microscope made in the U.S.A.



“Global microscopes are the foundation of my practice, as well as the reason for my comfort and efficiency.”

Dr. Stacey Ochoa
St. Louis, MO



“The Global microscope has changed my life and made my endodontics practice more enjoyable.”

Dr. Stewart Shapiro
Montreal, Quebec

CS | IMPLANTOLOGY | PERIODONTICS



2.1X



3.2X



5.1X

MAGNIFICATION CHART*

Magnification (x)/Field of View (dia-mm)

A3	3.2 / 62	6.4 / 32	12.8 / 16			
A4	3.2 / 62	5.1 / 39	8.0 / 25	12.8 / 16		
A6	2.1 / 95	3.2 / 62	5.1 / 39	8.0 / 25	12.8 / 16	19.2 / 10

*Using 10x eyepieces, 250mm objective lens, 160mm binocular.
With selected components, this magnification range can be shifted up or down.

VISUALIZATION | COMMUNICATION


EXCEPTIONAL COMMUNICATION

The A-Series™ microscope offers a phenomenal opportunity to improve communication between dental professionals, patients, and colleagues. Patients are given a chance to view live images, agree that there is a solution, and begin treatment options.

- Improved visualization through high levels of magnification
- Eyepieces produce exceptional field of view
- Diopter settings harmonize with eyeglasses
- Unobstructed, coaxial, shadow-free illumination
- Photos are exceptional means of documentation
- Deliver images with superior detail
- HD recording and display now available
- Use for presentations and teaching opportunities at conferences

“Microscopy is an integral part of how we diagnose our patients. Enlarged images and video are crucial for effective communication.”

Dr. Ronald Goldstein
Atlanta, GA



Photos courtesy of Dr. Glenn van As

N | PRESENTATION | CO-DIAGNOSIS



“After 20 years, I am still amazed at the impact the Global microscope has on my posture, precision, documentation, and communication.”

Dr. Glenn van As
Vancouver, BC

COMFORT | HEALTH | EFFICIENCY


LONG-TERM BENEFITS OF ERGONOMICS

Global microscopes enhance comfort and improve ergonomics for a more relaxed work day. Using the updated A-Series™ model promotes long-term health benefits by assisting the operator in maintaining optimum upright and balanced posture.

- Looking straight ahead through binoculars in a comfortable, stress-free position
- Ability to work longer periods without back and neck pain or fatigue
- Health benefits could extend practice life
- Increase productivity by seeing fewer patients while doing more dentistry on each one

“Within two weeks of using the Global microscope, my neck pain was gone. I have successfully managed my pain ever since.”

Dr. Lynn A. Jones
Bellevue, WA



LIFE-LONG HEALTH | Prevent back & neck strain

PRECISION | PRODUCTIVITY | POSTURE



GLOBAL

SURGICAL™ CORPORATION

We make it easier

800.767.8726 | globalsurgical.com

